

Going beyond  
the surface of  
your retina



**OCT-HS100**  
OPTICAL COHERENCE TOMOGRAPHY

*We Speak Image*

you can




**Canon**



**Automatic functions make examinations short and simple.**  
Perform the examination with only two simple mouse clicks!

#### 1. START



- Auto alignment
- Auto focus
- Auto C-gate

#### 2. CAPTURE



The Canon OCT-HS100 is the first Fully Automated Spectral Domain OCT. Operating an OCT has never been easier.

The OCT-HS100's extensive automatic functions simplify and optimize examinations. A single mouse click correctly aligns and focuses on the retina to maximize scan quality. The operator then decides with a second click when the scan itself will be initiated, allowing full control over the instrument. The high scan speed of 70,000 A-scans/sec results in very short examination times—typically under two seconds—improving efficiency and resulting in a very patient-friendly experience.

Canon's expertise in optics and innovative technology has resulted in a fantastic  $3\mu\text{m}$  optical axial resolution for amazing scan quality. The built in Scanning Laser Ophthalmoscope (SLO) allows for high quality retinal observations and precise follow up examinations.



# Full Auto OCT

High specifications in a very compact design

## Full Auto

The OCT-HS100's Full Auto feature significantly simplifies operation; standard procedures can easily be delegated and results are operator independent.

## High Image Quality.

3µm axial optical resolution

## Auto anterior eye alignment and tracking

The OCT-HS100 will automatically maintain the exact alignment on the center of the eye .

## Auto Fundus tracking by SLO

By detecting the amount of movement in fundus images, the unit can automatically compensate for small involuntary movements of the eye.

## Auto C-gate control (Coherence Gate)

Scan depth is automatically adjusted.

## Auto Focus

Automatic compensation of patient refraction.

## 10 mm scan width


Increased effectiveness of examinations by capturing large areas with just one scan.

## Easy examinations

Five stored examination sets: Macula Disease, Glaucoma, Choroid, Anterior and General. Examinations for other retinal diseases can be easily created in the initial settings menu.


## Extensive connectivity

The unit can be used stand-alone or in a network. Tomograms or reports can be output as JPEG, BMP or DICOM file.



## High scanning Speed


70,000 A-scans per second allow two-second examinations for reduced risk of motion artefacts and increased patient comfort.



## Built in SLO


For high quality retinal observation and precise follow up examinations

# Main Features




## 3 $\mu$ m axial resolution

3 micron resolution for unsurpassed image quality




## Averaging

Up to 50 scans can be combined for the highest possible image quality




## 10 layer recognition

The OCT-HS100 can recognize ten boundaries of the retina; detection of Bruch's membrane (BM) is now possible



## Choroid and Vitreous mode

Using Enhanced Depth Imaging, two modes can be selected based on the pathology: Choroid Mode (left) for detection of choroidal neovascularization. Vitreous Mode (below) for the detection of epiretinal membranes or vitreo macular traction;



# Extensive software



## Capture Screen

The capture screen displays a clear image of the anterior eye segment and retina (by SLO) as well as the OCT preview. The capture screen allows full control of the OCT-HS100: examination sets, fixation target, alignment and focusing.



## Macula analysis\*

Full retinal significance map  
Full retinal deviation map  
Full retinal difference map  
IS/OS - RPE thickness map  
RPE deformation map  
Inner retinal thickness map  
Outer retinal thickness map  
[NFL+GCL+IPL] thickness map

## ETDRS sectors


Full retinal measurements table  
Central sector thickness  
Minimum fovea thickness  
Total area average thickness  
Total volume  
Full retinal time-line graph  
3D object



## Macula glaucoma analysis\*

RNFL thickness map  
RNFL deviation map  
RNFL significance map  
ILM-RPE thickness map

NFL+GCL+IPL deviation map  
NFL+GCL+IPL significance map  
NFL+GCL+IPL sectors  
NFL+GCL+IPL time-line graph  
NFL+GCL+IPL measurement table



## Disc analysis\*

3.45mm circle reconstruction tomogram  
RNFL thickness  
RNFL deviation map  
RNFL significance map  
RNFL thickness profile  
ONH measurement table  
RNFL measurement table  
RNFL sectors  
RNFL time-line graph



## SLO tracking follow-up

SLO tracking allows accurate follow up examinations using the same scan position and identical scan parameters: Scan mode, position and size, fixation target position and size.

\* For these analyses a Normative Database will be available in a next version of the OCT-HS100 software

# Extensive scan modes



**Macula 3D**  
An area scan is done centering on the macula., for examination of the macula and the deeper layers  
Scan size 10 x 10 mm



**Cross**  
High resolution scan with up to 50 times averaging for highest image quality  
Scan size 3 x 3 to 10 x 10 mm



**Glaucoma 3D**  
A vertical area scan of the macula for extensive analysis of the NFL+GCL+IPL  
Scan size 10 x 10mm



**Multi Cross**  
Mult purpose high resolution scans with up to 10 times averaging.  
Scan size 3 x 3 to 10 x 10 mm




**Disc 3D**  
Area scan for disc analysis  
Scan size 6 x 6 mm



**Anterior 3D\***  
Area scan is performed on the anterior segment.  
Scan size 6 x 6 mm




**Custom**  
General application, for analysis of any disease  
Scan size 3 x 3 to 10 x 10 mm



**Anterior cross\***  
High resolution cross scan, for imaging cornea and angle.  
Scan size 3 x 3 to 10 x 10 mm

\* With optional adapter

## Retinal Image Import with automatic overlay



Imported retinal camera images are automatically aligned and overlaid on to the SLO image with precise position matching.

## Flexible Layout




The OCT-HS100 takes up very little floor space and is flexible for use in most situations—even against a wall or in a corner.

## Specifications


A-scans/sec	Max 70,000
Axial resolution	3 μm
Transversal Resolution	20 μm
Pupil size requirement	Min 3.0 mm
Scanning width	2 ~ 10 mm
Scan depth	2 mm
OCT light source	855 nm ±5 nm
Working Distance	35 mm

Fundus Preview	Fs SLO
Observation light source	780 ± 5nm
Internal Eye Fixation	2 mm or 6mm , 590nm (orange)
Field of view	10 x 10 mm, OCT 33° x 33°, SLO 44° x 33°
Dimensions (WxDxH)	387 x 499 x 474 (mm)
Weight	29 (kg)
Optional Accessory	Anterior segment adapter (ASA-1)


# Clinical images




Macular Edema



Macular Dystrophy



Posterior Epithelium Detachment



Age-related Macular Degeneration

Canon has been defining the future with innovative solutions for more than 70 years. In all that time we've constantly strived to improve medical diagnostics in healthcare. Perhaps that's what made us a leading global provider of eye care solutions.



#### Canon Eco

Our actions are based on honesty and sustainability.



#### Canon Quality

Safety and quality are an integral component of our actions.



#### Canon Flexibility

Everything we do has to have a superior customer advantage.

Choose the eye care system of the future and let our local, authorized Canon dealer advise you:

**Canon**

OCT-HS100  
English-NL Edition 2166V375  
© Canon Europa N.V. 2013

**Canon Europa N.V.  
Medical Systems Division**

Bovenkerkerweg 59-61  
1185 XB Amstelveen  
The Netherlands  
Phone: +31(0)20-545-8926  
Fax: +31(0)20-545-8220  
[www.canon-europe.com/medical](http://www.canon-europe.com/medical)